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Test strips with culture media to detect
aerobic bacteria, fungi and yeasts

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(Haltbarkeit 9 Monate).

Description

The CONTI+ hygiene tests are test strips for determining aerobic bacteria, fungi and yeasts in liquids and on surfaces. These test strips are a valuable aid for avoiding health risks from spray mist or droplet formation.

The hygiene tests allow easy testing of liquids and surfaces in order to determine the colony count (CFU). Within 48 hours the result regarding bacteria is available. A fungus detection can be made after 72 hours.

We recommend the use of an incubator for rapid cultivation / incubation of the medium at a temperature of about 30 °C.

The yellow surface of the hygiene test is intended for the detection of bacteria, the pink surface is for the detection of fungal infestation.

To read the result, compare the incubated hygiene test with the template.

Content

One packaging unit contains 10 hygiene test tubes with 2 immersion culture media each and note stickers for labelling the tubes.

Storage

The expiry date of the CONTI+ hygiene test is printed on the packaging.

The culture medium carriers are stored at room temperature (approx. 20 °C) protected from draughts and light, or in the refrigerator at 4 –8 °C. Condensation may occur with temperature fluctuations. However, this does not impair the quality of the CONTI+ Hygiene Test.

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
Usage

- Remove the culture medium carrier from the tube. Agar surfaces should not be touched.
- Immerse the agar surface in the liquid medium and wait for approx. 5 - 10 seconds.
- Drain excess liquid well and dab the lower edge of the support onto a clean filter paper or cloth.
- For surfaces, press each side of the culture medium only once on the surface to be tested. If possible, do not use yellow and pink sides on the same part of the surface.
- In the case of pasty or powder media, the agar surfaces should be coated with a sterile cotton swab. A separate swab is necessary for each agar surface.
- Afterwards the culture medium should be put back into the tube.
- There are note stickers in the packaging unit. Please label them (date, medium and place of collection) and then stick them on the tube.
- Incubation time of upright standing culture media: in incubator at approx. 30° C: for bacteria 48 hours, for yeasts and fungi 72 hours.
- At room temperature the incubation time is approx. 7 days at 20° C.

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Incubation

Incubation in warming cabinet
(27 °C – 30 °C)

for bacteria: 1–2 days
for yeasts and fungi: 3 days

Incubation at room temperature
(approx. 20 °C)

for bacteria: 2–3 days
for yeasts and fungi: 4–5 days

Analysis

Liquid samples are evaluated after the incubation period by comparing the colony density on the agar surfaces with the sample images (see figures).

The evaluation of contact samples is done by counting the colonies.

Notice:

If only very small colonies grow, it is recommended to extend the specified incubation period by 1–2 days.

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Evaluation:


up to 10^4 :	very weak to weakest contamination
from 10^4 - 10^5 :	moderate to severe contamination
above 10^6 :	strong to very strong contamination

Analysis

Determination of the total bacterial count (light yellow agar)

The addition of the indicator in the agar causes most bacterial germs to grow as red colonies. Colourless colonies must be taken into account when determining the number of bacteria. In rare cases, completely colourless growth or confluent bacterial growth may occur, which is shown as a uniformly colourless or red surface appears. In case of doubt, always compare the incubated immersion medium with an unused test strip. When assessing bacterial growth, the density of the colonies is relevant and not the size of the individual colonies.)

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Evaluation:

- + very weak to weakest contamination
- ++ moderate to severe contamination
- +++ strong to very strong contamination

Analysis

Determination of fungi (pink agar)

Filamentous fungi and yeasts grow on the pink agar. The growth can be pure thread fungus or yeast growth or mixed growth. Filamentous fungi form woolly colonies consisting of individual spores, filament parts or aggregates. Yeasts, on the other hand, form round, humpy and lustreless colonies, sometimes with runners.

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Evaluation:


up to 10^3 : very weak to weakest contamination

from $10^4 - 10^5$: moderate to severe contamination

above 10^6 : strong to very strong contamination

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Disposal

For proper disposal, we recommend placing contaminated immersion culture media overnight in a 0.15% disinfectant solution A-0803013. Alternatively, the culture media can also be incinerated or autoclaved.

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Reliable & innovative

For over 45 years, the CONTI+ brand has been a byword for flexible and holistic shower room and washroom solutions for public, semi-public and commercial environments as well as the health sector.

Known for its innovative technology and high grade materials, the brand has gained a reputation for unfailing quality. With a wide-ranging product portfolio and many years of experience creating custom designs, CONTI+ is now widely recognised among design engineers and architects for its valuable contribution to their projects. Every last detail is contemplated when developing tailor-made products, resulting in smart, individual solutions to suit a plethora of applications and equipment. The choice of shower and washbasin fittings encompasses a variety of colours and finishes as well as options for surface or flush mounting. Naturally, sustainability, hygiene and reliability are a prime consideration in the development of all products.

Thanks to its dependable quality and innovative technology, CONTI+ has remained a trusted brand in Germany for over 45 years.

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CONTI Sanitärarmaturen GmbH